

AS908 Rev.009

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

 **MB0926A**
**2026-11-10**

Instrument	Parameter		Low		Normal		High	
			 MB0926AL		 MB0926AN		 MB0926AH	
BC-7500[B] CRP	WBC	$\times 10^9/L$	3.76 $\pm$ 0.80		7.49 $\pm$ 1.00		20.81 $\pm$ 2.50	
BC-7500[N] CRP	Neu#	$\times 10^9/L$	2.13 $\pm$ 0.50		3.67 $\pm$ 0.80		11.28 $\pm$ 2.15	
BC-7500[R] CRP	Lym#	$\times 10^9/L$	0.91 $\pm$ 0.40		2.70 $\pm$ 0.80		6.80 $\pm$ 1.70	
BC-7500[NR] CRP	Mon#	$\times 10^9/L$	0.13 $\pm$ 0.20		0.37 $\pm$ 0.45		0.94 $\pm$ 1.15	
BC-7300[B] CRP	Eos#	$\times 10^9/L$	0.53 $\pm$ 0.40		0.65 $\pm$ 0.60		1.50 $\pm$ 1.50	
BC-7300[R] CRP	Bas#	$\times 10^9/L$	0.06 $\pm$ 0.22		0.10 $\pm$ 0.40		0.29 $\pm$ 1.00	
BC-7300[NR] CRP	Neu%	%	56.7 $\pm$ 12.0		48.9 $\pm$ 10.0		54.2 $\pm$ 10.0	
BC-7500[B] CS	Lym%	%	24.3 $\pm$ 9.5		36.0 $\pm$ 9.5		32.7 $\pm$ 8.0	
BC-7500[N] CS	Mon%	%	3.5 $\pm$ 5.0		5.0 $\pm$ 5.5		4.5 $\pm$ 5.5	
BC-7500[R] CS	Eos%	%	14.0 $\pm$ 9.0		8.7 $\pm$ 9.0		7.2 $\pm$ 8.0	
BC-7500[NR] CS	Bas%	%	1.5 $\pm$ 5.0		1.4 $\pm$ 5.0		1.4 $\pm$ 5.0	
QC Mode	RBC	$\times 10^{12}/L$	2.27 $\pm$ 0.18		4.27 $\pm$ 0.24		5.30 $\pm$ 0.30	
	HGB	g/dL	5.9 $\pm$ 0.4		12.3 $\pm$ 0.6		16.2 $\pm$ 0.8	
	HCT	%	19.1 $\pm$ 2.0		39.8 $\pm$ 2.5		53.0 $\pm$ 3.0	
	MCV	fL	84.2 $\pm$ 5.0		93.2 $\pm$ 5.0		100.0 $\pm$ 5.0	
	MCH	pg	26.0 $\pm$ 2.5		28.8 $\pm$ 2.5		30.6 $\pm$ 2.5	
	MCHC	g/dL	30.9 $\pm$ 3.0		30.9 $\pm$ 3.0		30.6 $\pm$ 3.0	
	RDW-CV	%	19.1 $\pm$ 3.0		15.3 $\pm$ 3.0		14.6 $\pm$ 3.0	
	RDW-SD	fL	58.7 $\pm$ 10.0		52.0 $\pm$ 10.0		53.1 $\pm$ 10.0	
	PLT	$\times 10^9/L$	54 $\pm$ 20		199 $\pm$ 45		405 $\pm$ 65	
	MPV	fL	9.7 $\pm$ 3.0		10.7 $\pm$ 3.0		11.0 $\pm$ 3.0	
	PDW	/	15.7 $\pm$ 3.0		16.6 $\pm$ 3.0		16.5 $\pm$ 3.0	
	PCT	%	0.052 $\pm$ 0.050		0.213 $\pm$ 0.100		0.446 $\pm$ 0.200	
	P-LCR	%	23.8 $\pm$ 20.0		31.4 $\pm$ 10.0		33.6 $\pm$ 10.0	
	P-LCC	$\times 10^9/L$	13 $\pm$ 12		62 $\pm$ 24		136 $\pm$ 50	
	IMG#	$\times 10^9/L$	0.05 $\pm$ 0.40		0.10 $\pm$ 0.80		0.27 $\pm$ 2.20	
	IMG%	%	1.3 $\pm$ 10.0		1.3 $\pm$ 10.0		1.3 $\pm$ 10.0	
	NRBC#	$\times 10^9/L$	0.318 $\pm$ 0.320		0.350 $\pm$ 0.320		0.489 $\pm$ 0.300	
	NRBC%	/100WBC	8.47 $\pm$ 8.25		4.67 $\pm$ 4.25		2.35 $\pm$ 2.30	

AS908 Rev.009

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

 **MB0926A**
**2026-11-10**

Instrument	Parameter		Low		Normal		High	
			 MB0926AL		 MB0926AN		 MB0926AH	
BC-6800 BC-6600 QC Mode	WBC	$\times 10^9/L$	3.77	± 0.80	7.31	± 1.00	20.25	± 2.50
	Neu#	$\times 10^9/L$	2.17	± 0.50	3.56	± 0.80	10.88	± 2.15
	Lym#	$\times 10^9/L$	0.90	± 0.40	2.67	± 0.80	6.72	± 1.70
	Mon#	$\times 10^9/L$	0.13	± 0.20	0.37	± 0.45	0.91	± 1.15
	Eos#	$\times 10^9/L$	0.53	± 0.40	0.63	± 0.60	1.52	± 1.50
	Bas#	$\times 10^9/L$	0.04	± 0.22	0.08	± 0.40	0.22	± 1.00
	Neu%	%	57.5	± 12.0	48.8	± 10.0	53.7	± 10.0
	Lym%	%	24.0	± 9.5	36.5	± 9.5	33.2	± 8.0
	Mon%	%	3.4	± 5.0	5.0	± 5.5	4.5	± 5.5
	Eos%	%	14.0	± 9.0	8.6	± 9.0	7.5	± 8.0
	Bas%	%	1.1	± 5.0	1.1	± 5.0	1.1	± 5.0
	RBC	$\times 10^{12}/L$	2.30	± 0.18	4.25	± 0.24	5.22	± 0.30
	HGB	g/dL	6.0	± 0.4	12.3	± 0.6	16.0	± 0.8
	HCT	%	19.6	± 2.0	39.7	± 2.5	52.0	± 3.0
	MCV	fL	85.1	± 5.0	93.4	± 5.0	99.6	± 5.0
	MCH	pg	26.1	± 2.5	28.9	± 2.5	30.7	± 2.5
	MCHC	g/dL	30.7	± 3.0	31.0	± 3.0	30.8	± 3.0
	RDW-CV	%	17.7	± 3.0	14.9	± 3.0	14.7	± 3.0
	RDW-SD	fL	53.4	± 10.0	49.8	± 10.0	51.7	± 10.0
	PLT	$\times 10^9/L$	53	± 20	199	± 45	416	± 65
	MPV	fL	9.4	± 3.0	10.3	± 3.0	10.7	± 3.0
	PDW	/	15.4	± 3.0	16.4	± 3.0	16.4	± 3.0
	PCT	%	0.050	± 0.050	0.205	± 0.100	0.445	± 0.200
	P-LCR	%	22.0	± 20.0	28.7	± 10.0	31.1	± 10.0
	P-LCC	$\times 10^9/L$	12	± 12	57	± 24	129	± 50
	IMG#*	$\times 10^9/L$	0.05	± 0.40	0.10	± 0.80	0.26	± 2.20
	IMG%*	%	1.3	± 10.0	1.3	± 10.0	1.3	± 10.0
	NRBC#	$\times 10^9/L$	0.375	± 0.320	0.394	± 0.320	0.563	± 0.300
	NRBC%	/100WBC	9.95	± 8.25	5.39	± 4.25	2.78	± 2.30

*Software version 1.09 or higher.

AS908 Rev.009

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

 **MB0926A**
**2026-11-10**

Instrument	Parameter		Low		Normal		High	
			LOT	MB0926AL	LOT	MB0926AN	LOT	MB0926AH
BC-6800Plus	WBC	×10 ⁹ /L	3.83	± 0.80	7.49	± 1.00	21.07	± 2.50
BC-6700Plus	Neu#	×10 ⁹ /L	2.19	± 0.50	3.65	± 0.80	11.31	± 2.15
BC-6600Plus	Lym#	×10 ⁹ /L	0.92	± 0.40	2.74	± 0.80	7.00	± 1.70
QC Mode	Mon#	×10 ⁹ /L	0.13	± 0.20	0.37	± 0.45	0.97	± 1.15
	Eos#	×10 ⁹ /L	0.53	± 0.40	0.63	± 0.60	1.50	± 1.50
	Bas#	×10 ⁹ /L	0.06	± 0.22	0.10	± 0.40	0.29	± 1.00
	Neu%	%	57.3	± 12.0	48.6	± 10.0	53.7	± 10.0
	Lym%	%	23.9	± 9.5	36.6	± 9.5	33.2	± 8.0
	Mon%	%	3.4	± 5.0	5.0	± 5.5	4.6	± 5.5
	Eos%	%	13.9	± 9.0	8.4	± 9.0	7.1	± 8.0
	Bas%	%	1.5	± 5.0	1.4	± 5.0	1.4	± 5.0
	RBC	×10 ¹² /L	2.34	± 0.18	4.29	± 0.24	5.26	± 0.30
	HGB	g/dL	6.0	± 0.4	12.2	± 0.6	15.9	± 0.8
	HCT	%	19.7	± 2.0	39.9	± 2.5	52.4	± 3.0
	MCV	fL	84.4	± 5.0	92.9	± 5.0	99.6	± 5.0
	MCH	pg	25.6	± 2.5	28.4	± 2.5	30.2	± 2.5
	MCHC	g/dL	30.4	± 3.0	30.6	± 3.0	30.3	± 3.0
	RDW-CV	%	17.8	± 3.0	15.1	± 3.0	14.7	± 3.0
	RDW-SD	fL	52.9	± 10.0	49.5	± 10.0	51.7	± 10.0
	PLT	×10 ⁹ /L	53	± 20	201	± 45	416	± 65
	PLT-H	×10 ⁹ /L	52	± 20	200	± 45	414	± 65
	MPV	fL	9.4	± 3.0	10.4	± 3.0	10.7	± 3.0
	PDW	/	15.4	± 3.0	16.3	± 3.0	16.3	± 3.0
	PCT	%	0.050	± 0.050	0.209	± 0.100	0.445	± 0.200
	P-LCR	%	21.5	± 20.0	29.3	± 10.0	31.5	± 10.0
	P-LCC	×10 ⁹ /L	11	± 12	59	± 24	131	± 50
	IMG#	×10 ⁹ /L	0.05	± 0.40	0.09	± 0.80	0.27	± 2.20
	IMG%	%	1.2	± 10.0	1.2	± 10.0	1.3	± 10.0
	NRBC#	×10 ⁹ /L	0.309	± 0.320	0.345	± 0.320	0.506	± 0.300
	NRBC%	/100WBC	8.08	± 8.25	4.60	± 4.25	2.40	± 2.30
Micro%	%	3.3	± 15.0	2.5	± 4.0	2.2	± 2.5	
Macro%	%	4.4	± 3.0	4.7	± 5.0	5.3	± 10.0	
Neu-X	/	379.4	± 120.0	378.5	± 120.0	379.8	± 120.0	
Neu-Y	/	505.4	± 120.0	508.8	± 150.0	510.3	± 150.0	
Mon-X	/	208.8	± 120.0	204.5	± 120.0	207.8	± 120.0	
PLT-I	×10 ⁹ /L	53	± 20	201	± 45	416	± 65	
ESR*	mm/h	37.97	± 10.00	16.09	± 4.00	/	± /	

*The ESR is applicable to software version 1.25 or higher, but not applicable for High level.

AS908 Rev.009

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

 **MB0926A**
**2026-11-10**

Instrument	Parameter		Low		Normal		High	
				MB0926AL		MB0926AN		MB0926AH
BC-6000	WBC	×10 ⁹ /L	3.65	± 0.80	7.38	± 1.00	20.88	± 2.50
BC-6100	Neu#	×10 ⁹ /L	2.10	± 0.50	3.58	± 0.80	11.18	± 2.15
BC-6200	Lym#	×10 ⁹ /L	0.87	± 0.40	2.69	± 0.80	6.95	± 1.70
Counter 60 AL	Mon#	×10 ⁹ /L	0.13	± 0.20	0.37	± 0.45	0.96	± 1.15
BC-6000Plus	Eos#	×10 ⁹ /L	0.50	± 0.40	0.63	± 0.60	1.48	± 1.50
BC-6100Plus	Bas#	×10 ⁹ /L	0.05	± 0.22	0.11	± 0.40	0.31	± 1.00
QC Mode	Neu%	%	57.3	± 12.0	48.6	± 10.0	53.5	± 10.0
	Lym%	%	23.9	± 9.5	36.4	± 9.5	33.3	± 8.0
	Mon%	%	3.5	± 5.0	5.0	± 5.5	4.6	± 5.5
	Eos%	%	13.8	± 9.0	8.5	± 9.0	7.1	± 8.0
	Bas%	%	1.5	± 5.0	1.5	± 5.0	1.5	± 5.0
	RBC	×10 ¹² /L	2.31	± 0.18	4.25	± 0.24	5.24	± 0.30
	HGB	g/dL	6.0	± 0.4	12.3	± 0.6	16.2	± 0.8
	HCT	%	20.0	± 2.0	40.0	± 2.5	52.1	± 3.0
	MCV	fL	86.6	± 5.0	94.0	± 5.0	99.4	± 5.0
	MCH	pg	26.0	± 2.5	28.9	± 2.5	30.9	± 2.5
	MCHC	g/dL	30.0	± 3.0	30.8	± 3.0	31.1	± 3.0
	RDW-CV	%	18.5	± 3.0	15.2	± 3.0	14.8	± 3.0
	RDW-SD	fL	57.7	± 10.0	52.1	± 10.0	53.4	± 10.0
	PLT	×10 ⁹ /L	54	± 20	207	± 45	427	± 65
	MPV	fL	10.5	± 3.0	11.1	± 3.0	11.2	± 3.0
	PDW	/	15.5	± 3.0	16.4	± 3.0	16.1	± 3.0
	PCT	%	0.057	± 0.050	0.230	± 0.100	0.478	± 0.200
	P-LCR	%	27.0	± 20.0	33.4	± 10.0	34.0	± 10.0
	P-LCC	×10 ⁹ /L	15	± 12	69	± 24	145	± 50
	IMG#	×10 ⁹ /L	0.05	± 0.40	0.10	± 0.80	0.27	± 2.20
	IMG%	%	1.3	± 10.0	1.3	± 10.0	1.3	± 10.0
	NRBC#	×10 ⁹ /L	0.315	± 0.320	0.348	± 0.320	0.484	± 0.300
	NRBC%	/100WBC	8.62	± 8.25	4.72	± 4.25	2.32	± 2.30
	Micro%	%	3.4	± 15.0	2.6	± 4.0	2.1	± 2.5
	Macro%	%	4.6	± 3.0	5.3	± 5.0	6.7	± 10.0
	Neu-X	/	398.9	± 120.0	386.9	± 120.0	397.2	± 120.0
	Neu-Y	/	518.6	± 120.0	518.6	± 150.0	516.7	± 150.0
	Mon-X	/	220.1	± 120.0	208.3	± 120.0	208.6	± 120.0
	PLT-I	×10 ⁹ /L	54	± 20	207	± 45	427	± 65

AS908 Rev.009

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT

 MB0926A


2026-11-10

Instrument	Parameter		Low		Normal		High	
			LOT	MB0926AL	LOT	MB0926AN	LOT	MB0926AH
BC-700[B]	WBC	$\times 10^9/L$		3.65 $\pm$ 0.80		7.52 $\pm$ 1.00		21.00 $\pm$ 2.50
BC-760[B]	Neu#	$\times 10^9/L$		2.17 $\pm$ 0.50		3.80 $\pm$ 0.80		11.69 $\pm$ 2.15
BC-700[B] CS	Lym#	$\times 10^9/L$		0.87 $\pm$ 0.40		2.73 $\pm$ 0.80		6.97 $\pm$ 1.70
BC-760[B] CS	Mon#	$\times 10^9/L$		0.13 $\pm$ 0.20		0.38 $\pm$ 0.45		0.97 $\pm$ 1.15
BC-700[R]	Eos#	$\times 10^9/L$		0.43 $\pm$ 0.40		0.50 $\pm$ 0.60		1.05 $\pm$ 1.50
BC-720[R]	Bas#	$\times 10^9/L$		0.05 $\pm$ 0.22		0.11 $\pm$ 0.40		0.32 $\pm$ 1.00
BC-760[R]	Neu%	%		59.4 $\pm$ 12.0		50.5 $\pm$ 10.0		55.7 $\pm$ 10.0
BC-700[R] CS	Lym%	%		23.7 $\pm$ 9.5		36.3 $\pm$ 9.5		33.2 $\pm$ 8.0
BC-760[R] CS	Mon%	%		3.5 $\pm$ 5.0		5.0 $\pm$ 5.5		4.6 $\pm$ 5.5
BP 200n	Eos%	%		11.9 $\pm$ 9.0		6.7 $\pm$ 9.0		5.0 $\pm$ 8.0
BP 260n	Bas%	%		1.5 $\pm$ 5.0		1.5 $\pm$ 5.0		1.5 $\pm$ 5.0
QC Mode	RBC	$\times 10^{12}/L$		2.29 $\pm$ 0.18		4.27 $\pm$ 0.24		5.26 $\pm$ 0.30
	HGB	g/dL		6.0 $\pm$ 0.4		12.3 $\pm$ 0.6		16.0 $\pm$ 0.8
	HCT	%		19.5 $\pm$ 2.0		40.3 $\pm$ 2.5		53.2 $\pm$ 3.0
	MCV	fL		85.0 $\pm$ 5.0		94.4 $\pm$ 5.0		101.2 $\pm$ 5.0
	MCH	pg		26.2 $\pm$ 2.5		28.8 $\pm$ 2.5		30.4 $\pm$ 2.5
	MCHC	g/dL		30.8 $\pm$ 3.0		30.5 $\pm$ 3.0		30.1 $\pm$ 3.0
	RDW-CV	%		19.0 $\pm$ 3.0		15.2 $\pm$ 3.0		14.3 $\pm$ 3.0
	RDW-SD	fL		59.3 $\pm$ 10.0		52.3 $\pm$ 10.0		53.2 $\pm$ 10.0
	PLT	$\times 10^9/L$		56 $\pm$ 20		203 $\pm$ 45		415 $\pm$ 65
	PLT-H	$\times 10^9/L$		57 $\pm$ 20		205 $\pm$ 45		413 $\pm$ 65
	MPV	fL		9.6 $\pm$ 3.0		10.5 $\pm$ 3.0		10.7 $\pm$ 3.0
	PDW	/		15.9 $\pm$ 3.0		16.6 $\pm$ 3.0		16.5 $\pm$ 3.0
	PCT	%		0.054 $\pm$ 0.050		0.213 $\pm$ 0.100		0.444 $\pm$ 0.200
	P-LCR	%		23.2 $\pm$ 20.0		30.6 $\pm$ 10.0		31.8 $\pm$ 10.0
	P-LCC	$\times 10^9/L$		13 $\pm$ 12		62 $\pm$ 24		132 $\pm$ 50
	IMG#	$\times 10^9/L$		0.05 $\pm$ 0.40		0.10 $\pm$ 0.80		0.27 $\pm$ 2.20
	IMG%	%		1.3 $\pm$ 10.0		1.3 $\pm$ 10.0		1.3 $\pm$ 10.0
	NRBC#	$\times 10^9/L$		0.085 $\pm$ 0.320		0.146 $\pm$ 0.320		0.456 $\pm$ 0.300
	NRBC%	/100WBC		2.33 $\pm$ 8.25		1.94 $\pm$ 4.25		2.17 $\pm$ 2.30
	IPF	%		2.4 $\pm$ 5.0		2.5 $\pm$ 5.0		2.4 $\pm$ 5.0
	ESR	mm/h		38.79 $\pm$ 10.00		14.85 $\pm$ 4.00		/ $\pm$ /

AS908 Rev.009

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT

 MB0926A


2026-11-10

Instrument	Parameter		Low		Normal		High	
			LOT	MB0926AL	LOT	MB0926AN	LOT	MB0926AH
BC-7600	WBC	$\times 10^9/L$	3.81	± 0.80	7.56	± 1.00	21.09	± 2.50
BC-7800	Neu#	$\times 10^9/L$	2.19	± 0.50	3.73	± 0.80	11.38	± 2.15
BC-7900 (Not applicable to Russia)	Lym#	$\times 10^9/L$	0.89	± 0.40	2.68	± 0.80	6.88	± 1.70
	Mon#	$\times 10^9/L$	0.13	± 0.20	0.37	± 0.45	0.97	± 1.15
	Eos#	$\times 10^9/L$	0.54	± 0.40	0.67	± 0.60	1.56	± 1.50
	Bas#	$\times 10^9/L$	0.06	± 0.22	0.11	± 0.40	0.30	± 1.00
	Neu%	%	57.4	± 12.0	49.4	± 10.0	54.0	± 10.0
	Lym%	%	23.4	± 9.5	35.4	± 9.5	32.6	± 8.0
	Mon%	%	3.4	± 5.0	4.9	± 5.5	4.6	± 5.5
	Eos%	%	14.3	± 9.0	8.9	± 9.0	7.4	± 8.0
	Bas%	%	1.5	± 5.0	1.4	± 5.0	1.4	± 5.0
	RBC	$\times 10^{12}/L$	2.27	± 0.18	4.24	± 0.24	5.24	± 0.30
	HGB	g/dL	6.0	± 0.4	12.3	± 0.6	16.1	± 0.8
	HCT	%	19.3	± 2.0	39.9	± 2.5	52.9	± 3.0
	MCV	fL	85.1	± 5.0	94.0	± 5.0	100.9	± 5.0
	MCH	pg	26.4	± 2.5	29.0	± 2.5	30.7	± 2.5
	MCHC	g/dL	31.1	± 3.0	30.9	± 3.0	30.5	± 3.0
	RDW-CV	%	19.3	± 3.0	15.4	± 3.0	14.6	± 3.0
	RDW-SD	fL	58.8	± 10.0	52.1	± 10.0	52.8	± 10.0
	PLT	$\times 10^9/L$	55	± 20	200	± 45	403	± 65
	PLT-H	$\times 10^9/L$	55	± 20	199	± 45	401	± 65
	MPV	fL	9.6	± 3.0	10.5	± 3.0	10.8	± 3.0
	PDW	/	15.8	± 3.0	16.6	± 3.0	16.5	± 3.0
	PCT	%	0.053	± 0.050	0.210	± 0.100	0.435	± 0.200
	P-LCR	%	23.4	± 20.0	30.5	± 10.0	32.6	± 10.0
	P-LCC	$\times 10^9/L$	13	± 12	61	± 24	131	± 50
	IMG#	$\times 10^9/L$	0.05	± 0.40	0.10	± 0.80	0.27	± 2.20
	IMG%	%	1.3	± 10.0	1.3	± 10.0	1.3	± 10.0
	NRBC#	$\times 10^9/L$	0.322	± 0.320	0.356	± 0.320	0.477	± 0.300
	NRBC%	/100WBC	8.45	± 8.25	4.71	± 4.25	2.26	± 2.30
	Micro%	%	3.1	± 15.0	2.4	± 4.0	2.0	± 2.5
	Macro%	%	4.5	± 3.0	4.8	± 5.0	5.6	± 10.0
	Neu-X	/	390.1	± 120.0	401.3	± 120.0	417.9	± 120.0
	Neu-Y	/	521.2	± 120.0	520.8	± 150.0	522.2	± 150.0
	Mon-X	/	224.3	± 120.0	215.7	± 120.0	220.3	± 120.0
	PLT-I	$\times 10^9/L$	55	± 20	200	± 45	403	± 65
	ESR	mm/h	38.36	± 10.00	14.19	± 4.00	/	$\pm /$
	WBC-M*	$\times 10^9/L$	3.82	± 0.80	7.53	± 1.00	21.35	± 2.50

The ESR is not applicable for High level, and only ESR is involved in BP 200n and BP 260n.

WBC-M applies to BC-7900.

*Research use only.

AS908 Rev.009