

AS908 Rev.010

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT

 MB1126A


2027-01-10

Instrument	Parameter		Low		Normal		High	
			LOT	MB1126AL	LOT	MB1126AN	LOT	MB1126AH
BC-7500[B] CRP	WBC	$\times 10^9/L$		3.70 $\pm$ 0.80		7.65 $\pm$ 1.00		20.81 $\pm$ 2.50
BC-7500[N] CRP	Neu#	$\times 10^9/L$		2.06 $\pm$ 0.50		3.75 $\pm$ 0.80		11.27 $\pm$ 2.15
BC-7500[R] CRP	Lym#	$\times 10^9/L$		0.91 $\pm$ 0.40		2.64 $\pm$ 0.80		6.85 $\pm$ 1.70
BC-7500[NR] CRP	Mon#	$\times 10^9/L$		0.13 $\pm$ 0.20		0.37 $\pm$ 0.45		0.96 $\pm$ 1.15
BC-7300[B] CRP	Eos#	$\times 10^9/L$		0.54 $\pm$ 0.40		0.78 $\pm$ 0.60		1.42 $\pm$ 1.50
BC-7300[R] CRP	Bas#	$\times 10^9/L$		0.06 $\pm$ 0.22		0.11 $\pm$ 0.40		0.31 $\pm$ 1.00
BC-7300[NR] CRP	Neu%	%		55.7 $\pm$ 12.0		48.9 $\pm$ 10.0		54.2 $\pm$ 10.0
BC-7500[B] CS	Lym%	%		24.6 $\pm$ 9.5		34.5 $\pm$ 9.5		32.9 $\pm$ 8.0
BC-7500[N] CS	Mon%	%		3.5 $\pm$ 5.0		4.9 $\pm$ 5.5		4.6 $\pm$ 5.5
BC-7500[R] CS	Eos%	%		14.7 $\pm$ 9.0		10.2 $\pm$ 9.0		6.8 $\pm$ 8.0
BC-7500[NR] CS	Bas%	%		1.5 $\pm$ 5.0		1.5 $\pm$ 5.0		1.5 $\pm$ 5.0
QC Mode	RBC	$\times 10^{12}/L$		2.26 $\pm$ 0.18		4.20 $\pm$ 0.24		5.16 $\pm$ 0.30
	HGB	g/dL		5.8 $\pm$ 0.4		12.2 $\pm$ 0.6		15.8 $\pm$ 0.8
	HCT	%		18.9 $\pm$ 2.0		39.4 $\pm$ 2.5		51.9 $\pm$ 3.0
	MCV	fL		83.5 $\pm$ 5.0		93.8 $\pm$ 5.0		100.6 $\pm$ 5.0
	MCH	pg		25.7 $\pm$ 2.5		29.0 $\pm$ 2.5		30.6 $\pm$ 2.5
	MCHC	g/dL		30.7 $\pm$ 3.0		31.0 $\pm$ 3.0		30.4 $\pm$ 3.0
	RDW-CV	%		19.0 $\pm$ 3.0		15.4 $\pm$ 3.0		14.5 $\pm$ 3.0
	RDW-SD	fL		57.5 $\pm$ 10.0		52.6 $\pm$ 10.0		52.3 $\pm$ 10.0
	PLT	$\times 10^9/L$		54 $\pm$ 20		198 $\pm$ 45		412 $\pm$ 65
	MPV	fL		11.1 $\pm$ 3.0		10.7 $\pm$ 3.0		9.5 $\pm$ 3.0
	PDW	/		16.9 $\pm$ 3.0		16.5 $\pm$ 3.0		16.6 $\pm$ 3.0
	PCT	%		0.060 $\pm$ 0.050		0.212 $\pm$ 0.100		0.391 $\pm$ 0.200
	P-LCR	%		37.1 $\pm$ 20.0		30.7 $\pm$ 10.0		23.6 $\pm$ 10.0
	P-LCC	$\times 10^9/L$		20 $\pm$ 12		61 $\pm$ 24		97 $\pm$ 50
	IMG#	$\times 10^9/L$		0.05 $\pm$ 0.40		0.10 $\pm$ 0.80		0.27 $\pm$ 2.20
	IMG%	%		1.3 $\pm$ 10.0		1.3 $\pm$ 10.0		1.3 $\pm$ 10.0
	NRBC#	$\times 10^9/L$		0.312 $\pm$ 0.320		0.354 $\pm$ 0.320		0.425 $\pm$ 0.300
	NRBC%	/100WBC		8.43 $\pm$ 8.25		4.63 $\pm$ 4.25		2.04 $\pm$ 2.30

AS908 Rev.010

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

 **MB1126A**
**2027-01-10**

Instrument	Parameter		Low		Normal		High	
			 MB1126AL		 MB1126AN		 MB1126AH	
BC-6800 BC-6600 QC Mode	WBC	$\times 10^9/L$	3.81	± 0.80	7.53	± 1.00	20.32	± 2.50
	Neu#	$\times 10^9/L$	2.12	± 0.50	3.67	± 0.80	11.03	± 2.15
	Lym#	$\times 10^9/L$	0.96	± 0.40	2.61	± 0.80	6.64	± 1.70
	Mon#	$\times 10^9/L$	0.13	± 0.20	0.37	± 0.45	0.91	± 1.15
	Eos#	$\times 10^9/L$	0.56	± 0.40	0.80	± 0.60	1.52	± 1.50
	Bas#	$\times 10^9/L$	0.04	± 0.22	0.08	± 0.40	0.22	± 1.00
	Neu%	%	55.4	± 12.0	48.8	± 10.0	54.2	± 10.0
	Lym%	%	25.3	± 9.5	34.7	± 9.5	32.7	± 8.0
	Mon%	%	3.5	± 5.0	4.9	± 5.5	4.5	± 5.5
	Eos%	%	14.7	± 9.0	10.6	± 9.0	7.5	± 8.0
	Bas%	%	1.1	± 5.0	1.0	± 5.0	1.1	± 5.0
	RBC	$\times 10^{12}/L$	2.30	± 0.18	4.22	± 0.24	5.14	± 0.30
	HGB	g/dL	5.9	± 0.4	12.3	± 0.6	15.8	± 0.8
	HCT	%	19.3	± 2.0	39.6	± 2.5	51.4	± 3.0
	MCV	fL	83.8	± 5.0	93.8	± 5.0	100.0	± 5.0
	MCH	pg	25.7	± 2.5	29.1	± 2.5	30.7	± 2.5
	MCHC	g/dL	30.6	± 3.0	31.1	± 3.0	30.7	± 3.0
	RDW-CV	%	17.7	± 3.0	15.1	± 3.0	14.6	± 3.0
	RDW-SD	fL	52.2	± 10.0	50.2	± 10.0	51.8	± 10.0
	PLT	$\times 10^9/L$	53	± 20	200	± 45	412	± 65
	MPV	fL	11.0	± 3.0	10.4	± 3.0	9.4	± 3.0
	PDW	/	16.7	± 3.0	16.1	± 3.0	16.5	± 3.0
	PCT	%	0.058	± 0.050	0.208	± 0.100	0.387	± 0.200
	P-LCR	%	36.1	± 20.0	28.5	± 10.0	22.6	± 10.0
	P-LCC	$\times 10^9/L$	19	± 12	57	± 24	93	± 50
	IMG#*	$\times 10^9/L$	0.05	± 0.40	0.10	± 0.80	0.26	± 2.20
	IMG%*	%	1.3	± 10.0	1.3	± 10.0	1.3	± 10.0
	NRBC#	$\times 10^9/L$	0.350	± 0.320	0.378	± 0.320	0.492	± 0.300
	NRBC%	/100WBC	9.19	± 8.25	5.02	± 4.25	2.42	± 2.30

*Software version 1.09 or higher.

AS908 Rev.010

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

 **MB1126A**
**2027-01-10**

Instrument	Parameter	Low		Normal		High	
		 MB1126AL		 MB1126AN		 MB1126AH	
BC-6800Plus	WBC	$\times 10^9/L$	3.91 $\pm$ 0.80	7.82 $\pm$ 1.00		21.44 $\pm$ 2.50	
BC-6700Plus	Neu#	$\times 10^9/L$	2.16 $\pm$ 0.50	3.75 $\pm$ 0.80		11.45 $\pm$ 2.15	
BC-6600Plus	Lym#	$\times 10^9/L$	0.98 $\pm$ 0.40	2.78 $\pm$ 0.80		7.16 $\pm$ 1.70	
QC Mode	Mon#	$\times 10^9/L$	0.14 $\pm$ 0.20	0.38 $\pm$ 0.45		0.99 $\pm$ 1.15	
	Eos#	$\times 10^9/L$	0.57 $\pm$ 0.40	0.80 $\pm$ 0.60		1.54 $\pm$ 1.50	
	Bas#	$\times 10^9/L$	0.06 $\pm$ 0.22	0.11 $\pm$ 0.40		0.30 $\pm$ 1.00	
	Neu%	%	55.4 $\pm$ 12.0	47.9 $\pm$ 10.0		53.4 $\pm$ 10.0	
	Lym%	%	25.1 $\pm$ 9.5	35.6 $\pm$ 9.5		33.4 $\pm$ 8.0	
	Mon%	%	3.5 $\pm$ 5.0	4.9 $\pm$ 5.5		4.6 $\pm$ 5.5	
	Eos%	%	14.5 $\pm$ 9.0	10.2 $\pm$ 9.0		7.2 $\pm$ 8.0	
	Bas%	%	1.5 $\pm$ 5.0	1.4 $\pm$ 5.0		1.4 $\pm$ 5.0	
	RBC	$\times 10^{12}/L$	2.35 $\pm$ 0.18	4.29 $\pm$ 0.24		5.17 $\pm$ 0.30	
	HGB	g/dL	5.9 $\pm$ 0.4	12.3 $\pm$ 0.6		15.7 $\pm$ 0.8	
	HCT	%	19.6 $\pm$ 2.0	40.2 $\pm$ 2.5		52.0 $\pm$ 3.0	
	MCV	fL	83.3 $\pm$ 5.0	93.7 $\pm$ 5.0		100.5 $\pm$ 5.0	
	MCH	pg	25.1 $\pm$ 2.5	28.7 $\pm$ 2.5		30.4 $\pm$ 2.5	
	MCHC	g/dL	30.1 $\pm$ 3.0	30.6 $\pm$ 3.0		30.2 $\pm$ 3.0	
	RDW-CV	%	17.8 $\pm$ 3.0	15.1 $\pm$ 3.0		14.4 $\pm$ 3.0	
	RDW-SD	fL	52.3 $\pm$ 10.0	50.1 $\pm$ 10.0		51.0 $\pm$ 10.0	
	PLT	$\times 10^9/L$	53 $\pm$ 20	199 $\pm$ 45		410 $\pm$ 65	
	PLT-H	$\times 10^9/L$	52 $\pm$ 20	198 $\pm$ 45		408 $\pm$ 65	
	MPV	fL	10.9 $\pm$ 3.0	10.4 $\pm$ 3.0		9.4 $\pm$ 3.0	
	PDW	/	16.5 $\pm$ 3.0	16.1 $\pm$ 3.0		16.5 $\pm$ 3.0	
	PCT	%	0.058 $\pm$ 0.050	0.207 $\pm$ 0.100		0.385 $\pm$ 0.200	
	P-LCR	%	35.8 $\pm$ 20.0	28.8 $\pm$ 10.0		22.6 $\pm$ 10.0	
	P-LCC	$\times 10^9/L$	19 $\pm$ 12	57 $\pm$ 24		93 $\pm$ 50	
	IMG#	$\times 10^9/L$	0.05 $\pm$ 0.40	0.10 $\pm$ 0.80		0.26 $\pm$ 2.20	
	IMG%	%	1.2 $\pm$ 10.0	1.3 $\pm$ 10.0		1.2 $\pm$ 10.0	
	NRBC#	$\times 10^9/L$	0.329 $\pm$ 0.320	0.363 $\pm$ 0.320		0.474 $\pm$ 0.300	
	NRBC%	/100WBC	8.42 $\pm$ 8.25	4.64 $\pm$ 4.25		2.21 $\pm$ 2.30	
	Micro%	%	3.7 $\pm$ 15.0	2.5 $\pm$ 4.0		2.1 $\pm$ 2.5	
	Macro%	%	4.1 $\pm$ 3.0	4.9 $\pm$ 5.0		5.4 $\pm$ 10.0	
	Neu-X	/	398.5 $\pm$ 120.0	393.3 $\pm$ 120.0		401.7 $\pm$ 120.0	
	Neu-Y	/	519.4 $\pm$ 120.0	517.8 $\pm$ 150.0		519.2 $\pm$ 150.0	
	Mon-X	/	216.5 $\pm$ 120.0	225.1 $\pm$ 120.0		227.1 $\pm$ 120.0	
	PLT-I	$\times 10^9/L$	53 $\pm$ 20	199 $\pm$ 45		410 $\pm$ 65	
	ESR*	mm/h	36.69 $\pm$ 10.00	15.87 $\pm$ 4.00		/ $\pm$ /	

*The ESR is applicable to software version 1.25 or higher, but not applicable for High level.

AS908 Rev.010

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

 **MB1126A**
**2027-01-10**

Instrument	Parameter		Low		Normal		High	
				MB1126AL		MB1126AN		MB1126AH
BC-6000	WBC	×10 ⁹ /L	3.65	± 0.80	7.56	± 1.00	20.61	± 2.50
BC-6100	Neu#	×10 ⁹ /L	2.04	± 0.50	3.66	± 0.80	11.17	± 2.15
BC-6200	Lym#	×10 ⁹ /L	0.90	± 0.40	2.62	± 0.80	6.78	± 1.70
Counter 60 AL	Mon#	×10 ⁹ /L	0.13	± 0.20	0.37	± 0.45	0.95	± 1.15
BC-6000Plus	Eos#	×10 ⁹ /L	0.53	± 0.40	0.80	± 0.60	1.40	± 1.50
BC-6100Plus	Bas#	×10 ⁹ /L	0.05	± 0.22	0.11	± 0.40	0.31	± 1.00
QC Mode	Neu%	%	55.8	± 12.0	48.3	± 10.0	54.2	± 10.0
	Lym%	%	24.7	± 9.5	34.7	± 9.5	32.9	± 8.0
	Mon%	%	3.5	± 5.0	4.9	± 5.5	4.6	± 5.5
	Eos%	%	14.5	± 9.0	10.6	± 9.0	6.8	± 8.0
	Bas%	%	1.5	± 5.0	1.5	± 5.0	1.5	± 5.0
	RBC	×10 ¹² /L	2.31	± 0.18	4.26	± 0.24	5.15	± 0.30
	HGB	g/dL	6.0	± 0.4	12.4	± 0.6	15.8	± 0.8
	HCT	%	19.8	± 2.0	40.0	± 2.5	51.2	± 3.0
	MCV	fL	85.8	± 5.0	94.0	± 5.0	99.4	± 5.0
	MCH	pg	26.0	± 2.5	29.1	± 2.5	30.7	± 2.5
	MCHC	g/dL	30.3	± 3.0	31.0	± 3.0	30.9	± 3.0
	RDW-CV	%	18.2	± 3.0	15.3	± 3.0	14.3	± 3.0
	RDW-SD	fL	56.5	± 10.0	53.3	± 10.0	53.3	± 10.0
	PLT	×10 ⁹ /L	55	± 20	205	± 45	434	± 65
	MPV	fL	12.1	± 3.0	11.3	± 3.0	9.7	± 3.0
	PDW	/	16.6	± 3.0	16.1	± 3.0	16.5	± 3.0
	PCT	%	0.067	± 0.050	0.232	± 0.100	0.421	± 0.200
	P-LCR	%	42.2	± 20.0	33.5	± 10.0	24.0	± 10.0
	P-LCC	×10 ⁹ /L	23	± 12	69	± 24	104	± 50
	IMG#	×10 ⁹ /L	0.05	± 0.40	0.10	± 0.80	0.27	± 2.20
	IMG%	%	1.3	± 10.0	1.3	± 10.0	1.3	± 10.0
	NRBC#	×10 ⁹ /L	0.311	± 0.320	0.355	± 0.320	0.410	± 0.300
	NRBC%	/100WBC	8.51	± 8.25	4.70	± 4.25	1.99	± 2.30
	Micro%	%	3.6	± 15.0	2.5	± 4.0	2.0	± 2.5
	Macro%	%	4.6	± 3.0	5.5	± 5.0	6.8	± 10.0
	Neu-X	/	399.6	± 120.0	400.2	± 120.0	405.1	± 120.0
	Neu-Y	/	517.1	± 120.0	516.9	± 150.0	521.5	± 150.0
	Mon-X	/	217.6	± 120.0	209.3	± 120.0	209.7	± 120.0
	PLT-I	×10 ⁹ /L	55	± 20	205	± 45	434	± 65

AS908 Rev.010

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

 **MB1126A**
**2027-01-10**

Instrument	Parameter		Low		Normal		High	
			 MB1126AL		 MB1126AN		 MB1126AH	
BC-700[B]	WBC	$\times 10^9/L$	3.74	± 0.80	7.84	± 1.00	21.33	± 2.50
BC-760[B]	Neu#	$\times 10^9/L$	2.15	± 0.50	3.96	± 0.80	12.03	± 2.15
BC-700[B] CS	Lym#	$\times 10^9/L$	0.92	± 0.40	2.72	± 0.80	6.97	± 1.70
BC-760[B] CS	Mon#	$\times 10^9/L$	0.13	± 0.20	0.38	± 0.45	0.96	± 1.15
BC-700[R]	Eos#	$\times 10^9/L$	0.48	± 0.40	0.66	± 0.60	1.05	± 1.50
BC-720[R]	Bas#	$\times 10^9/L$	0.06	± 0.22	0.12	± 0.40	0.32	± 1.00
BC-760[R]	Neu%	%	57.5	± 12.0	50.5	± 10.0	56.4	± 10.0
BC-700[R] CS	Lym%	%	24.7	± 9.5	34.7	± 9.5	32.7	± 8.0
BC-760[R] CS	Mon%	%	3.5	± 5.0	4.9	± 5.5	4.5	± 5.5
BP 200n	Eos%	%	12.8	± 9.0	8.4	± 9.0	4.9	± 8.0
BP 260n	Bas%	%	1.5	± 5.0	1.5	± 5.0	1.5	± 5.0
QC Mode	RBC	$\times 10^{12}/L$	2.29	± 0.18	4.24	± 0.24	5.15	± 0.30
	HGB	g/dL	5.9	± 0.4	12.3	± 0.6	15.8	± 0.8
	HCT	%	19.1	± 2.0	40.0	± 2.5	52.2	± 3.0
	MCV	fL	83.5	± 5.0	94.3	± 5.0	101.3	± 5.0
	MCH	pg	25.8	± 2.5	29.0	± 2.5	30.7	± 2.5
	MCHC	g/dL	30.9	± 3.0	30.8	± 3.0	30.3	± 3.0
	RDW-CV	%	18.8	± 3.0	15.2	± 3.0	14.2	± 3.0
	RDW-SD	fL	57.9	± 10.0	52.8	± 10.0	52.7	± 10.0
	PLT	$\times 10^9/L$	56	± 20	205	± 45	414	± 65
	PLT-H	$\times 10^9/L$	56	± 20	202	± 45	408	± 65
	MPV	fL	11.2	± 3.0	10.7	± 3.0	9.3	± 3.0
	PDW	/	17.0	± 3.0	16.5	± 3.0	16.6	± 3.0
	PCT	%	0.063	± 0.050	0.219	± 0.100	0.385	± 0.200
	P-LCR	%	36.6	± 20.0	29.9	± 10.0	22.7	± 10.0
	P-LCC	$\times 10^9/L$	20	± 12	61	± 24	94	± 50
	IMG#	$\times 10^9/L$	0.05	± 0.40	0.10	± 0.80	0.28	± 2.20
	IMG%	%	1.3	± 10.0	1.3	± 10.0	1.3	± 10.0
	NRBC#	$\times 10^9/L$	0.084	± 0.320	0.151	± 0.320	0.471	± 0.300
	NRBC%	/100WBC	2.24	± 8.25	1.93	± 4.25	2.21	± 2.30
	IPF	%	2.6	± 5.0	2.4	± 5.0	2.4	± 5.0
	ESR	mm/h	37.61	± 10.00	14.63	± 4.00	/	$\pm /$

AS908 Rev.010

BC-6D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

 **MB1126A**
**2027-01-10**

Instrument	Parameter		Low		Normal		High	
			LOT	MB1126AL	LOT	MB1126AN	LOT	MB1126AH
BC-7600	WBC	×10 ⁹ /L	3.85	± 0.80	7.81	± 1.00	21.27	± 2.50
BC-7800	Neu#	×10 ⁹ /L	2.15	± 0.50	3.83	± 0.80	11.68	± 2.15
BC-7900	Lym#	×10 ⁹ /L	0.96	± 0.40	2.72	± 0.80	6.93	± 1.70
BC-7600Plus	Mon#	×10 ⁹ /L	0.13	± 0.20	0.37	± 0.45	0.98	± 1.15
BC-7800Plus	Eos#	×10 ⁹ /L	0.55	± 0.40	0.78	± 0.60	1.38	± 1.50
BC-7900Plus	Bas#	×10 ⁹ /L	0.06	± 0.22	0.11	± 0.40	0.30	± 1.00
(Not applicable to Russia)	Neu%	%	55.7	± 12.0	49.0	± 10.0	54.9	± 10.0
	Lym%	%	24.9	± 9.5	34.8	± 9.5	32.6	± 8.0
	Mon%	%	3.5	± 5.0	4.8	± 5.5	4.6	± 5.5
	Eos%	%	14.4	± 9.0	10.0	± 9.0	6.5	± 8.0
	Bas%	%	1.5	± 5.0	1.4	± 5.0	1.4	± 5.0
	RBC	×10 ¹² /L	2.31	± 0.18	4.26	± 0.24	5.17	± 0.30
	HGB	g/dL	5.9	± 0.4	12.3	± 0.6	15.8	± 0.8
	HCT	%	19.4	± 2.0	40.3	± 2.5	52.4	± 3.0
	MCV	fL	83.9	± 5.0	94.5	± 5.0	101.3	± 5.0
	MCH	pg	25.5	± 2.5	28.9	± 2.5	30.6	± 2.5
	MCHC	g/dL	30.4	± 3.0	30.6	± 3.0	30.2	± 3.0
	RDW-CV	%	19.3	± 3.0	15.6	± 3.0	14.6	± 3.0
	RDW-SD	fL	57.8	± 10.0	52.9	± 10.0	52.7	± 10.0
	PLT	×10 ⁹ /L	58	± 20	203	± 45	416	± 65
	PLT-H	×10 ⁹ /L	57	± 20	202	± 45	414	± 65
	MPV	fL	10.9	± 3.0	10.6	± 3.0	9.3	± 3.0
	PDW	/	16.9	± 3.0	16.4	± 3.0	16.7	± 3.0
	PCT	%	0.063	± 0.050	0.215	± 0.100	0.387	± 0.200
	P-LCR	%	36.2	± 20.0	29.9	± 10.0	23.0	± 10.0
	P-LCC	×10 ⁹ /L	21	± 12	61	± 24	96	± 50
	IMG#	×10 ⁹ /L	0.05	± 0.40	0.10	± 0.80	0.28	± 2.20
	IMG%	%	1.3	± 10.0	1.3	± 10.0	1.3	± 10.0
	NRBC#	×10 ⁹ /L	0.326	± 0.320	0.364	± 0.320	0.436	± 0.300
	NRBC%	/100WBC	8.46	± 8.25	4.66	± 4.25	2.05	± 2.30
	Micro%	%	3.4	± 15.0	2.4	± 4.0	2.0	± 2.5
	Macro%	%	4.5	± 3.0	5.0	± 5.0	5.6	± 10.0
	Neu-X	/	408.1	± 120.0	406.7	± 120.0	405.7	± 120.0
	Neu-Y	/	522.4	± 120.0	522.2	± 150.0	525.6	± 150.0
	Mon-X	/	221.0	± 120.0	225.4	± 120.0	223.1	± 120.0
	PLT-I	×10 ⁹ /L	58	± 20	203	± 45	416	± 65
	ESR	mm/h	37.33	± 10.00	14.43	± 4.00	/	± /
	WBC-M*	×10 ⁹ /L	3.83	± 0.80	7.96	± 1.00	21.66	± 2.50

The ESR is not applicable for High level, and only ESR is involved in BP 200n and BP 260n.

WBC-M applies to BC-7900/BC-7900Plus.

*Research use only.

AS908 Rev.010